

Attention: Users of this C6 Flow Cytometer®

Before running samples:

- 1) Place an empty 12x75 mm tube on the SIP.
- 2) Click **BACKFLUSH**.
- 3) Place a **fresh** tube with 2 mL of filtered, deionized H₂O on the SIP.
- 4) Set time limit for 2 minutes and Fluidics speed to **Fast**.
- 5) Click **RUN**.
- 6) Once time limit is reached, click **DELETE SAMPLE DATA**.
- 7) Remove the tube and run samples.

After running samples:

- 1) Place a tube with 2 mL of diluted Cleaning Solution (#KR-210) on the SIP*.
- 2) Select an empty data well.
- 3) Set time limit for 2 minutes and Fluidics speed to **Fast**.
- 4) Click **RUN**.
- 5) Once time limit is reached, remove the tube of Cleaning Solution.
- 6) Place a tube of 2 mL filtered, DI H₂O on the SIP.
- 7) Set a time limit of two 2 minutes.
- 8) Click **RUN**. The C6 will stop automatically when the time limit is reached.
- 9) **Leave the tube on the SIP until the C6 is used again.**

*If a more aggressive cleaning is needed, Decontamination Solution (KR-200) may be used in place of Cleaning Solution.

DO NOT USE BD FACSTM Rinse or BD FACSTM Clean or other non-Accuri approved cleaning or decontamination fluids.

7820023A